



PRODUCT DESCRIPTION

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS OR ANIMALS

ProteinTech Group, Inc
2201 W. Campbell Park Dr.
Chicago, IL 60612
Toll Free: 1-888-4PTGLAB
In Illinois: 1312-455-8498
Fax: 1-312-455-8408
Web: www.ptglab.com


Antibody to GFP

The antibody:

Format 1: Affinity purified
Catalog No. : 50430-2-AP
Weight/Volume: 300ug/ml
Tested Applications: WB; IHC; IP; IF

The production

Immunogen: recombinant protein
Tested Specificity: Hu, Mo
Host: Rb
Antibody species: Polyclonal



293 cell lysates transfected with pGFP vector of fusion protein annexin A5 were subjected to SDS PAGE followed by western blot with affinity purified Ab at dilution of 1:1000-3000

Storage Information:

The antibody is stored at 50% Glycerol in PBS buffer containing 0.02% sodium azide at -20 °C. The final concentration generally is 150 µg/ml for affinity purified antibody. Glycerol was added to a final concentration of 50% so that antibodies can be stored at -20 °C without going through freeze-thaw cycles.