



PRODUCT DESCRIPTION

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS OR ANIMALS

ProteinTech Group, Inc
2201 W. Campbell Park Dr.
Chicago, IL 60612
Toll Free: 1-888-4PTGLAB
In Illinois: 1312-455-8498
Fax: 1-312-455-8408
Web: www.ptglab.com

Antibody to Flag

The antibody:


Format 1: Affinity purified
Catalog No. : 50144-3-AP
Weight/Volume: 400 µg/ml
Tested Applications: WB; IHC; IP; IF

The production

Immunogen: peptide
Tested Specificity: Hu, Mo
Host: Rb
Antibody species: Polyclonal

Storage Information:

The antibody is stored at 50% Glycerol in PBS buffer containing 0.02% sodium azide at -20 °C. The final concentration generally is 150 µg/ml for affinity purified antibody. Glycerol was added to a final concentration of 50% so that antibodies can be stored at -20 °C without going through freeze-thaw cycles.



293 cell lysates transfected with p3XFLAG-CMV vector of fusion protein TBC1D7 were subjected to SDS PAGE followed by western blot with affinity purified Ab at dilution of 1:1000-3000