

For Research Use Only

APC Armenian Hamster IgG Isotype Control (PIP)



Catalog Number: APC-65210

2 Publications

Basic Information

Catalog Number:

APC-65210

Size:

100ug, 0.2 mg/mL

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Affinity purification

CloneNo.:

PIP

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Background Information

The PIP immunoglobulin is useful as an isotype-matched control. The PIP immunoglobulin is specific for Glutathione S-transferase and can be used as an isotype control for Armenian Hamster IgG antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Notable Publications

Author	Pubmed ID	Journal	Application
Qian Gou	37870169	Adv Sci (Weinh)	FC
Qian Gou	37096255	Cancer Sci	

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.25 ug APC Armenian Hamster IgG Isotype Control (APC-65210, Clone: P1P) (blue) or 0.25 ug APC Anti-Mouse CD3 ϵ (APC-65060, Clone: 145-2C11) (red). Cells were not fixed.