

For Research Use Only

APC Rat IgG1 Isotype Control (HRPN)

Catalog Number: APC-65212



Basic Information

Catalog Number:

APC-65212

Concentration:

100ug, 0.2 mg/mL

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Affinity purification

CloneNo.:

HRPN

Recommended Dilutions:

FC: 0.5 ug per 10⁶ cells in 100 μ L suspension

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

Positive Controls:

FC :,

Background Information

The HRPN immunoglobulin is useful as an isotype-matched control. The HRPN immunoglobulin has an unknown binding specificity and is used as an isotype control for rat IgG1 antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.1% gelatin, pH7.2

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.5 ug APC Rat IgG1 Isotype Control (APC-65212, Clone: HRPN) (blue), or not stained (blank) (black). Cells were not fixed.