

内参与与标签抗体



CoraLite®-Conjugated

PCNA HRP- β -Actin

GAPDH Lamin B1 BrdU

ATPIA1 6*his, His tag COX4I2

α -Tubulin β -Tubulin COX4I1

GFP tag DYKDDDDK tag

HA Tag FITC-conjugated

Proteintech

为人类的每一个蛋白开发出相应的抗体



旗下品牌

Proteintech已完成针对人类基因靶标的抗体 SCI引用文献超过17万次

· 内参抗体

适用范围	名称	分子量	抗体类型	货号	应用	WB 推荐稀释比	文献
胞浆或全细胞	GAPDH	36 kDa	鼠单抗	60004-1-Ig	WB, IF, FC, IP, IHC, ChIP, CoIP	1:5000-1:50000	5979
			兔多抗	10494-1-AP	IHC, WB, IP, IF, FC, CoIP, RIP	1:5000-1:40000	3442
	β-Actin	42 kDa	鼠单抗	66009-1-Ig	WB, IF, FC, IP, IHC, ChIP, CoIP	1:2000-1:100000	3106
			兔多抗	20536-1-AP	FC, IF, IHC, WB, IP	1:1000-1:5000	1650
	α-Tubulin	50-55 kDa	鼠单抗	66031-1-Ig	IHC, FC, IP, WB, IF	1:2000-1:100000	482
			兔多抗	11224-1-AP	WB, IP, IHC, FC, IF, CoIP	1:2000-1:16000	755
细胞核膜			重组兔单抗	80762-1-RR	IF, WB	1:5000-1:50000	—
			鼠单抗	66240-1-Ig	WB, IHC, IF, IP	1:2000-1:100000	166
	β-Tubulin	50-55 kDa	兔多抗	10094-1-AP	WB, IF, IP, IHC, FC, CoIP	1:1000-1:5000	415
			重组兔单抗	80713-1-RR	FC, IF, IP, WB	1:5000-1:50000	—
	Vinculin	117 kDa	鼠单抗	66305-1-Ig	WB, IF, IP, IHC, FC, CoIP	1:5000-1:10000	59
	Lamin A/C	65-75 kDa	兔多抗	10298-1-AP	IP, FC, WB, IF	1:1000-1:6000	125
细胞核	Lamin B1	66-72 kDa	鼠单抗	66095-1-Ig	WB, FC, IP, IF, IHC	1:1000-1:50000	191
			兔多抗	12987-1-AP	IHC, IP, WB, FC, IF, ChIP	1:2000-1:10000	601
	TBP	33-43 kDa*	鼠单抗	66166-1-Ig	WB, IP, IHC, IF, FC	1:2000-1:16000	16
			兔多抗	22006-1-AP	WB, IP, IF	1:500-1:3000	50
	PCNA	36 kDa	鼠单抗	60097-1-Ig	IHC, FC, IP, WB, IF	1:5000-1:50000	58
			兔多抗	10205-2-AP	IF, IP, WB, FC, IHC	1:2000-1:10000	525
膜蛋白	Histone-H3	15-17 kDa	鼠单抗	17168-1-AP	WB, IHC, IF, IP, FC, ChIP	1:2000-1:16000	490
	YY1	65-70 kDa	鼠单抗	66281-1-Ig	WB, IF, IP, IHC, ChIP, CoIP	1:5000-1:50000	18
			兔多抗	22156-1-AP	ChIP, IF, IHC, IP, WB	1:500-1:2000	19
	ATP1A1	97-110 kDa	兔多抗	14418-1-AP	IF, IHC, IP, WB	1:5000-1:20000	100
			兔多抗	55187-1-AP	WB, IHC, FC, IF	1:1000-1:8000	9
			鼠单抗	66345-1-Ig	WB, IHC, IF	1:2000-1:20000	9
线粒体	VDAC1/Porin	31-37 kDa	兔多抗	10866-1-AP	CoIP, IF, IHC, WB	1:500-1:3000	98
	COX411	17-20 kDa	鼠单抗	66110-1-Ig	FC, IF, IHC, WB	1:2000-1:8000	14
			兔多抗	11242-1-AP	FC, IF, IHC, IP, WB, ChIP	1:5000-1:50000	255
	COX412	17-20 kDa	兔多抗	11463-1-AP	WB, IHC, IF, IP	1:500-1:3000	15
	BrdU	—	鼠单抗	66241-1-Ig	IHC, IF	—	47
	Transferrin	77 kDa	兔多抗	66171-1-Ig	IHC, WB, IF	1:5000-1:50000	2
全血 / 血浆 / 血清			兔多抗	17435-1-AP	WB, IHC, IF	1:1000-1:5000	24
	Albumin	66 kDa	鼠单抗	66051-1-Ig	WB, IP, IHC, IF	1:5000-1:50000	15
			兔多抗	16475-1-AP	WB, IP, IHC, IF, CoIP	1:5000-1:50000	75

*TBP 在人类是 37-43 kDa, 在小鼠和大鼠中是 33-36 kDa

· 标签抗体

名称	类型	货号	应用	WB 推荐稀释比	文献
6*his,His tag	兔多抗	10001-0-AP	IF, IP, WB	1:500-1:2000	23
	鼠单抗	66005-1-Ig	ChIP, CoIP, ELISA, FC, IF, IHC, IP, WB	1:5000-1:50000	447
DYKDDDDK tag	鼠单抗	66008-4-Ig	ChIP, CoIP, IP, WB, IF	1:5000-1:50000	6
	兔多抗	20543-1-AP	ChIP, CoIP, FC, IF, IHC, IP, RIP, WB	1:2000-1:16000	431
GFP tag	鼠单抗	66002-1-Ig	ChIP, CoIP, IF, IHC, IP, WB, FC	1:2000-1:16000	266
	兔多抗	50430-2-AP	CoIP, IF, IHC, IP, IP-MS, RIP, WB	1:1000-1:4000	350
GST Tag	鼠单抗	66001-2-Ig	CoIP, IF, IP, WB	1:5000-1:50000	40
	兔多抗	10000-0-AP	ChIP, IF, IP, WB, FC	1:1000-1:4000	113
HA Tag	鼠单抗	66006-2-Ig	CoIP, FC, IF, IP, PLA, WB	1:10000-1:100000	79
	兔多抗	51064-2-AP	ChIP, CoIP, IF, IHC, IP, WB	1:5000-1:10000	339
MYC tag	鼠单抗	60003-2-Ig	CoIP, IF, IP, WB	1:2000-1:10000	87
	兔多抗	16286-1-AP	CoIP, FC, IF, IHC, IP, WB	1:1000-1:8000	127
V5-tag	兔多抗	14440-1-AP	IF, IP, WB	1:2000-1:20000	15
S1 tag	鼠单抗	66165-1-Ig	WB	1:5000-1:50000	—
MBP Tag	鼠单抗	66003-1-Ig	CoIP, IF, IHC, IP, WB	1:1000-1:8000	22
	兔多抗	15089-1-AP	IF, IHC, WB	1:1000-1:6000	14

· 标记内参 / 标签抗体

类别	鼠单抗名称	货号	应用	WB 推荐稀释比	文献
HRP 标记内参抗体	HRP-anti GAPDH	HRP-60004	IF, IHC, WB	1:5000-1:50000	341
	HRP-anti β-Actin	HRP-66009	IP, WB	1:2000-1:50000	23
	HRP-anti α-Tubulin	HRP-66031	WB	1:2000-1:16000	62
HRP 标记标签抗体	HRP-anti 6*His, His-Tag	HRP-66005	Blocking Assay, WB	1:5000-1:50000	70
	HRP-anti GFP tag	HRP-66002	IF, WB	1:5000-1:50000	6
	HRP-anti GST Tag	HRP-66001	IP, WB	1:5000-1:50000	19
	HRP-conjugated DYKDDDDK Tag	HRP-66008	WB	1:10000-1:100000	7
	CoraLite™488-Conjugated Beta Actin	CL488-66009	FC, IF, WB	WB : 1:2000-1:10000	—
	CoraLite™594-Conjugated Beta Actin	CL594-66009	FC (Intra), IF, WB	WB : 1:500-1:2000	2
荧光标记内参抗体	CoraLite™532-Conjugated Beta Actin	CL532-66009	IF	—	—
	CoraLite™555-Conjugated Beta Actin	CL555-66009	FC (Intra), IF	—	—
	CoraLite™647-Conjugated Beta Actin	CL647-66009	FC	—	—
	CoraLite™488-Conjugated Alpha Tubulin	CL488-66031	IF	—	—
	CoraLite™594-Conjugated Alpha Tubulin	CL594-66031	IF	—	—
	CoraLite™555-Conjugated Alpha Tubulin	CL555-66031	FC (Intra), IF	—	—
	CoraLite™647-Conjugated Alpha Tubulin	CL647-66031	FC (Intra)	—	—
	CoraLite™488-conjugated 6*His, His-Tag	CL488-66005	IF	—	1
	CoraLite™594-conjugated 6*His, His-Tag	CL594-66005	IF	—	—
	CoraLite™647-conjugated 6*His, His-Tag	CL647-66005	IF	—	—
荧光标记标签抗体	FITC-conjugated DYKDDDDK tag	FITC-66008	IF	—	—
	CoraLite™488-conjugated DYKDDDDK tag	CL488-80010	IF	—	—
	CoraLite™594-conjugated DYKDDDDK tag	CL594-80010	IF	—	—
	CoraLite™647-conjugated DYKDDDDK tag	CL647-80010	IF	—	—
	CoraLite™488-conjugated GFP tag	CL488-66002	WB	1:2000-1:10000	—
	CoraLite™594-conjugated GFP tag	CL594-66002	FC, IF, WB	1:2000-1:10000	—
	CoraLite™647-conjugated GFP tag	CL647-66002	IF	—	—
	CoraLite™488-conjugated GST tag	CL488-66001	WB	1:5000-1:50000	—
	CoraLite™594-conjugated GST tag	CL594-50430	IF, WB	—	—
	CoraLite™647-conjugated GST tag	CL647-10000	FC	—	—
	CoraLite™488-conjugated HA tag	CL488-66003	IF	—	—
	CoraLite™594-conjugated HA tag	CL594-51064	IF	—	—
	CoraLite™647-conjugated HA tag	CL647-60003	IF	—	—