

For Research Use Only

# K2-15mer Polyclonal antibody

Catalog Number: 10040-1-AP



---

## Basic Information

Catalog Number:  
10040-1-AP  
Size:  
N/A  
Source:  
Rabbit  
Isotype:  
IgG

GenBank Accession Number:  
GeneID (NCBI):  
Full Name:  
Calculated MW:  
0 kDa

Purification Method:  
Antigen affinity purification

---

## Applications

Tested Applications:  
ELISA  
Species Specificity:  
human

---

## Background Information

---

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

---

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data