

For Research Use Only

SIN3A Polyclonal antibody

Catalog Number: 10046-1-AP **1 Publications**



Basic Information

Catalog Number:

10046-1-AP

Concentration:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066364

GeneID (NCBI):

25942

UNIPROT ID:

Q96ST3

Full Name:

SIN3 homolog A, transcription regulator (yeast)

Calculated MW:

145 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse

Cited Species:

human, pig

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Mingyang Li	39826728	Int J Biol Macromol	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data