

For Research Use Only

Calmodulin1/2/3 Polyclonal antibody



Catalog Number: 10112-1-AP

Basic Information

Catalog Number:

10112-1-AP

Size:

507 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0146

GenBank Accession Number:

BC005137

GeneID (NCBI):

808

UNIPROT ID:

P62158

Full Name:

calmodulin 3 (phosphorylase kinase, delta)

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human lung cancer tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung cancer using 10112-1-AP (CALM3 antibody) at dilution of 1:50 (under 10x lens).