

For Research Use Only

THRA Polyclonal antibody

Catalog Number: 10139-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number:

10139-1-AP

Size:

2500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0182

GenBank Accession Number:

BC002728

GeneID (NCBI):

7067

UNIPROT ID:

P10827

Full Name:

thyroid hormone receptor, alpha
(erythroblastic leukemia viral (v-erb-
a) oncogene homolog, avian)

Calculated MW:

490 aa, 54 kDa

Observed MW:

58 kDa, 50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

Positive Controls:

WB: A431 cells, mouse brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sijie Wang	36216077	Physiol Behav	WB
C B Steinhauser	34062290	Domest Anim Endocrinol	IHC
Jun Guo	39639130	Hypertens Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 10139-1-AP (THRA antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.