

For Research Use Only

# DAP1 Polyclonal antibody

Catalog Number: 10160-1-AP



## Basic Information

Catalog Number:

10160-1-AP

Size:

287 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0181

GenBank Accession Number:

BC002726

GeneID (NCBI):

1611

UNIPROT ID:

P51397

Full Name:

death-associated protein

Calculated MW:

16 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IF, IHC, ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 10160-1-AP (DAP1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using 10160-1-AP (DAP1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).