

For Research Use Only

TAP2 Polyclonal antibody, PBS Only

Catalog Number:10161-1-PBS



Basic Information

Catalog Number:

10161-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0190

GenBank Accession Number:

BC002751

GeneID (NCBI):

6891

UNIPROT ID:

Q03519

Full Name:

transporter 2, ATP-binding cassette, sub-family B (MDR/TAP)

Calculated MW:

76 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

TAP2, also known as APT2, ABCB3, is a member of the superfamily of ATP-binding cassette (ABC) transporters. TAP2 mediates unidirectional translocation of peptide antigens from cytosol to endoplasmic reticulum (ER) for loading onto MHC class I (MHCI) molecules (PMID: 25656091). Mutations in TAP2 may be associated with Bare lymphocyte syndrome 1 (BLS1) and schizophrenia (PMID: 7517574, 24365204).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10161-1-AP (TAP2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 10161-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human intrahepatic cholangiocarcinoma tissue slide using 10161-1-AP (TAP2 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 10161-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 10161-1-AP (TAP2 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 10161-1-PBS in a different storage buffer formulation.