

For Research Use Only

MAGEA3 Polyclonal antibody

Catalog Number: 10220-1-AP

1 Publications



Basic Information

Catalog Number:

10220-1-AP

Size:

100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0265

GenBank Accession Number:

BC000340

GeneID (NCBI):

4102

UNIPROT ID:

P43357

Full Name:

melanoma antigen family A, 3

Calculated MW:

314 aa, 35 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Kangchan Choi	39410046	Cancers (Basel)	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

