

For Research Use Only

# PRPF18 Polyclonal antibody

Catalog Number: 10261-1-AP



## Basic Information

Catalog Number:

10261-1-AP

Concentration:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0329

GenBank Accession Number:

BC002572

GeneID (NCBI):

8559

UNIPROT ID:

Q99633

Full Name:

PRP18 pre-mRNA processing factor 18  
homolog (S. cerevisiae)

Calculated MW:

5 aa, 475 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, L02 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10261-1-AP (PRPF18 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.