

For Research Use Only

# PIR Polyclonal antibody

Catalog Number:10263-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

10263-1-AP

Size:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0330

GenBank Accession Number:

BC002517

GeneID (NCBI):

8544

UNIPROT ID:

O00625

Full Name:

pirin (iron-binding nuclear protein)

Calculated MW:

32 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

## Applications

Tested Applications:

WB, IHC, IP, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HepG2 cells, HeLa cells

IP : HeLa cells,

IHC : human stomach tissue,

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

## Background Information

Pirin (iron-binding nuclear protein, PIR) is a highly conserved 32-kDa protein consisting of 290 amino acids. Pirin mRNA is expressed weakly in all human tissues. Confocal immunofluorescence experiments demonstrate a predominant localization of Pirin within dot-like subnuclear structures. Pirin interacts with the oncoprotein B-cell lymphoma 3-encoded (Bcl-3) and nuclear factor I (NFI).

## Notable Publications

Author	Pubmed ID	Journal	Application
Muhammad Suleman	31500513	Cell Cycle	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 10263-1-AP (PIR antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



IP result of anti-PIR (IP:10263-1-AP, 4ug; Detection:10263-1-AP 1:300) with HeLa cells lysate 1200ug.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 10263-1-AP (PIR antibody) at dilution of 1:100 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).