

For Research Use Only

CCNI Polyclonal antibody

Catalog Number: 10279-1-AP **1 Publications**



Basic Information

Catalog Number: 10279-1-AP	GenBank Accession Number: BC000420	Purification Method: Antigen affinity purification
Size: 547 µg/ml	GeneID (NCBI): 10983	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	UNIPROT ID: Q14094	
Isotype: IgG	Full Name: cyclin I	
Immunogen Catalog Number: AG0300	Calculated MW: 43 kDa	

Applications

Tested Applications: IF/ICC,ELISA	Positive Controls: IF : HepG2 cells,
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Wentao Zhu	38356330	J Med Chem	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data