

For Research Use Only

# Goat-IgG Polyclonal antibody

Catalog Number: 10285-1-AP

50 Publications



## Basic Information

Catalog Number:

10285-1-AP

Size:

1400 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

155 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF, RIP

Species Specificity:

goat

Cited Species:

mouse

## Background Information

### Notable Publications

| Author       | Pubmed ID | Journal                   | Application |
|--------------|-----------|---------------------------|-------------|
| Jian Wen     | 32996777  | Cancer Biother Radiopharm | WB          |
| Lan Zhang    | 34586738  | Clin Transl Med           | RIP         |
| Hucheng Chen | 30272275  | Mol Med Rep               |             |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

