

For Research Use Only

Calmodulin 1/2/3 Polyclonal antibody



Catalog Number: 10303-1-AP

Basic Information

Catalog Number:

10303-1-AP

Size:

100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0490

GenBank Accession Number:

BC000454

GeneID (NCBI):

801

UNIPROT ID:

PODP23

Full Name:

calmodulin 1 (phosphorylase kinase, delta)

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung cancer using 10303-1-AP (CALM1 antibody) at dilution of 1:100 (under 25x lens).