

For Research Use Only

SMA4 Polyclonal antibody

Catalog Number: 10318-1-AP



Basic Information

Catalog Number: 10318-1-AP	GenBank Accession Number: BC002622	Purification Method: Antigen affinity purification
Size: 600 µg/ml	GeneID (NCBI): 11039	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: glucuronidase, beta pseudogene	
Isotype: IgG	Calculated MW: 15 kDa	
Immunogen Catalog Number: AG0372	Observed MW: 15 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 10318-1-AP (SMA4 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.