

For Research Use Only

PPP2R4 Polyclonal antibody

Catalog Number: 10321-1-AP

Featured Product

9 Publications



Basic Information

Catalog Number:

10321-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0392

GenBank Accession Number:

BC002545

GeneID (NCBI):

5524

UNIPROT ID:

Q15257

Full Name:

protein phosphatase 2A activator,
regulatory subunit 4

Calculated MW:

37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:200-1:800

Applications

Tested Applications:

IHC, WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB : HEK-293 cells, HeLa cells, HepG2 cells, Jurkat
cells, K-562 cells

IHC : mouse testis tissue, mouse kidney tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Liping Chen	25236601	J Cell Sci	
Mingming Yang	35648484	Nucleic Acids Res	WB
Xinxing Pan	35545004	Mol Cancer Ther	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of PP2A antibody (10321-1-AP, 1:600) with si-control and si-PP2A transfected HEK293 cells.



Various lysates were subjected to SDS PAGE followed by western blot with 10321-1-AP (PPP2R4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 10321-1-AP (PPP2R4 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 10321-1-AP (PPP2R4 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).