

For Research Use Only

# Melanophilin Polyclonal antibody

Catalog Number: 10338-1-AP

Featured Product

19 Publications



## Basic Information

### Catalog Number:

10338-1-AP

### Size:

800 µg/ml

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG0316

### GenBank Accession Number:

BC001653

### GeneID (NCBI):

79083

### UNIPROT ID:

Q9BV36

### Full Name:

melanophilin

### Calculated MW:

66 kDa

### Observed MW:

66 kDa, 80 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

IF 1:20-1:200

## Applications

### Tested Applications:

IF/ICC, IHC, IP, WB, ELISA

### Cited Applications:

WB, IP, IF, IHC

### Species Specificity:

human

### Cited Species:

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

### Positive Controls:

WB : A375 cells,

IP : A375 cells,

IHC : human skin cancer tissue,

IF : HeLa cells, MCF-7 cells

## Background Information

Melanophilin (MLph) is a member of the exophilin subfamily of Rab effector proteins. MLPH is a modular protein, forming a ternary complex with the small Ras-related GTPase Rab27A in its GTP-bound form and the motor protein myosin Va. It regulates retention of melanosomes at the peripheral actin cytoskeleton of melanocytes, a process essential for normal mammalian pigmentation. Defects in MLPH are a cause of Griscelli syndrome type 3 (GS3), which is characterized by a silver-gray hair color and abnormal pigment distribution in the hair shaft.

## Notable Publications

Author	Pubmed ID	Journal	Application
Valentina Cetica	25312756	J Allergy Clin Immunol	WB, IP
Jacob Wang	19468668	J Gastroenterol	IHC
Qiaorong Huang	35602484	Front Med (Lausanne)	WB

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data




A375 cells were subjected to SDS PAGE followed by western blot with 10338-1-AP (Melanophilin antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 10338-1-AP (Melanophilin antibody at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 10338-1-AP (Melanophilin antibody at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using 10338-1-AP (Melanophilin antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



IP result of anti-Melanophilin (IP:10338-1-AP, 4ug; Detection:10338-1-AP 1:600) with A375 cells lysate 3600ug.