

For Research Use Only

NGFRAP1 Polyclonal antibody

Catalog Number: 10446-1-AP



Basic Information

Catalog Number:

10446-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0378

GenBank Accession Number:

BC003190

GeneID (NCBI):

27018

UNIPROT ID:

Q00994

Full Name:

nerve growth factor receptor
(TNFRSF16) associated protein 1

Calculated MW:

13 kDa, 22 kDa, 44 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HepG2 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 lysates were subjected to SDS PAGE followed by western blot with 10446-1-AP (NGFRAP1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.