

For Research Use Only

ZMAT3 Polyclonal antibody

Catalog Number: 10504-1-AP

8 Publications



Basic Information

Catalog Number:

10504-1-AP

Concentration:

1100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0780

GenBank Accession Number:

BC002896

GeneID (NCBI):

64393

UNIPROT ID:

Q9HA38

Full Name:

zinc finger, matrin type 3

Calculated MW:

32 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

WB, IP, ELISA

Cited Applications:

WB, IHC, IF, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Positive Controls:

WB : COLO 320 cells, rat testis tissue, HEK-293 cells, HepG2 cells, mouse testis tissue

IP : mouse testis tissue,

Background Information

ZMAT3 is a p53-induced zinc finger protein, and a role in the TP53-dependent growth regulatory pathway, contribute to TP53-mediated apoptosis by regulation of TP53 expression and translocation to the nucleus and nucleolus. Acts as a bona fide target gene of p53/TP53.

Notable Publications

Author	Pubmed ID	Journal	Application
Yuan Zhang	30833572	Cell Death Dis	WB
Xin Yang	35850611	J Physiol Sci	WB
Jian Zhang	33439503	J Cell Mol Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10504-1-AP (ZMAT3 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 10504-1-AP (ZMAT3 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.



IP result of anti-ZMAT3 (IP:10504-1-AP, 3ug; Detection:10504-1-AP 1:700) with mouse testis tissue lysate 4000ug.