

For Research Use Only

# NIF3L1 Polyclonal antibody

Catalog Number: 10549-1-AP



## Basic Information

Catalog Number:

10549-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0819

GenBank Accession Number:

BC007654

GeneID (NCBI):

60491

UNIPROT ID:

Q9GZT8

Full Name:

NIF3 NGG1 interacting factor 3-like 1  
(S. pombe)

Calculated MW:

42 kDa

Observed MW:

42 kDa, 39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : HeLa cells, HEK-293 cells, HepG2 cells, PC-3 cells,  
PC-12 cells, C6 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10549-1-AP (NIF3L1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.