

For Research Use Only

NIF3L1 Polyclonal antibody

Catalog Number: 10549-1-AP



Basic Information

Catalog Number:

10549-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0819

GenBank Accession Number:

BC007654

GeneID (NCBI):

60491

UNIPROT ID:

Q9GZT8

Full Name:

NIF3 NGG1 interacting factor 3-like 1
(S. pombe)

Calculated MW:

42 kDa

Observed MW:

42 kDa, 39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : HeLa cells, HEK-293 cells, HepG2 cells, PC-3 cells,
PC-12 cells, C6 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10549-1-AP (NIF3L1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.