

For Research Use Only

# BL21-lysate Polyclonal antibody

Catalog Number: 10557-1-Ig



## Basic Information

Catalog Number:

10557-1-Ig

Size:

1800  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0049

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

1-300 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:1000-1:8000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

bacteria

Positive Controls:

WB: Total bacteria lysates of BL21,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

