

For Research Use Only

BL21(DE3)-lysate Polyclonal antibody



Catalog Number: 10558-1-Ig **1 Publications**

Basic Information

Catalog Number:

10558-1-Ig

Size:

2000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0048

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

1-300 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

bacteria

Cited Species:

zebrafish

Positive Controls:

WB: Total bacteria lysates of BL21(DE3),

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sakura Tanaka	34978328	Endocrinology	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data