

For Research Use Only

ARL6IP4 Polyclonal antibody

Catalog Number: 10563-1-AP



Basic Information

Catalog Number:

10563-1-AP

Size:

587 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0863

GenBank Accession Number:

BC001958

GeneID (NCBI):

51329

UNIPROT ID:

Q66PJ3

Full Name:

ADP-ribosylation-like factor 6
interacting protein 4

Calculated MW:

25 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

ARL6IP4 is recognized to be a binding protein of ARL6 that is another member in the Ras family, implicating a possibility of mediating Ras signaling. ARL6 gene mutation has been linked to Bardet-Biedl syndrome type 3 which is a multisystemic disorder characterized by defects in the microtubule-based transport process.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

