

For Research Use Only

SPOCK3 Polyclonal antibody

Catalog Number: 10603-1-AP



Basic Information

Catalog Number:

10603-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0838

GenBank Accession Number:

BC003017

GeneID (NCBI):

50859

UNIPROT ID:

Q9BQ16

Full Name:

sparc/osteonectin, cwcv and kazal-like domains proteoglycan (testican)
3

Calculated MW:

49 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

