

For Research Use Only

ERC1 Polyclonal antibody

Catalog Number: 10619-1-AP



Basic Information

Catalog Number: 10619-1-AP	GenBank Accession Number: BC005065	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 23085	
Source: Rabbit	UNIPROT ID: Q8IUD2	
Isotype: IgG	Full Name: ELKS/RAB6-interacting/CAST family member 1	
Immunogen Catalog Number: AG0966	Calculated MW: 128 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data