

For Research Use Only

IFRD2 Polyclonal antibody

Catalog Number: 10649-1-AP



Basic Information

Catalog Number: 10649-1-AP	GenBank Accession Number: BC007265	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 7866	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q12894	
Isotype: IgG	Full Name: interferon-related developmental regulator 2	
Immunogen Catalog Number: AG1021	Calculated MW: 48 kDa	
	Observed MW: 50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, A549 cells, PC-3 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 10649-1-AP (IFRD2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.