

For Research Use Only

ACTRT3/ARPM1 Polyclonal antibody

Catalog Number: 10652-1-AP



Basic Information

Catalog Number:

10652-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1010

GenBank Accession Number:

BC007289

GeneID (NCBI):

84517

UNIPROT ID:

Q9BYD9

Full Name:

actin related protein M1

Calculated MW:

41 kDa

Observed MW:

40-43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HT-29 cells, U-87 MG cells, mouse brain tissue, mouse kidney tissue, mouse testis tissue, rat brain tissue, rat kidney tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10652-1-AP (ARPM1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.