

For Research Use Only

UBL5 Polyclonal antibody, PBS Only

Catalog Number:10715-1-PBS



Basic Information

Catalog Number:

10715-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1122

GenBank Accession Number:

BC007053

GeneID (NCBI):

59286

UNIPROT ID:

Q9BZL1

Full Name:

ubiquitin-like 5

Calculated MW:

8.5 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse heart tissue slide using 10715-1-AP (UBL5 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 10715-1-PBS in a different storage buffer formulation.