

For Research Use Only

TCEAL3 Polyclonal antibody

Catalog Number: 10732-2-AP



Basic Information

Catalog Number: 10732-2-AP	GenBank Accession Number: BC007622	Purification Method: Antigen affinity purification
Size: 150 µg/ml	GeneID (NCBI): 85012	
Source: Rabbit	UNIPROT ID: Q969E4	
Isotype: IgG	Full Name: transcription elongation factor A (SII)-like 3	
Immunogen Catalog Number: AG1039	Calculated MW: 23 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human ovary tumor using 10732-2-AP (TCEAL3 antibody) at dilution of 1:100 (under 40x lens).