

For Research Use Only

# TCEAL3 Polyclonal antibody

Catalog Number: 10732-2-AP



## Basic Information

Catalog Number:

10732-2-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1039

GenBank Accession Number:

BC007622

GeneID (NCBI):

85012

UNIPROT ID:

Q969E4

Full Name:

transcription elongation factor A (SII)-like 3

Calculated MW:

23 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human ovary tumor using 10732-2-AP (TCEAL3 antibody) at dilution of 1:100 (under 40x lens).