

For Research Use Only

Aurora B Polyclonal antibody

Catalog Number: 10792-1-AP **1 Publications**



Basic Information

Catalog Number: 10792-1-AP	GenBank Accession Number: BC013300	Purification Method: Antigen affinity purification
Size: 133 µg/ml	GeneID (NCBI): 9212	
Source: Rabbit	UNIPROT ID: Q96GD4	
Isotype: IgG	Full Name: aurora kinase B	
Immunogen Catalog Number: AG1246	Calculated MW: 39.2 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Aurora B(AURKB) is also named as aurora kinase B, AIK2, AIM1, AIRK2, ARK2, STK1, STK12, STK5 and belongs to the Aurora subfamily. Aurora kinase B is a chromosomal passenger serine/threonine protein kinase that regulates accurate chromosomal segregation, cytokinesis, protein localization to the centromere and kinetochore, correct microtubule-kinetochore attachments, and regulation of the mitotic checkpoint(PMID:18311747). Expression levels of AURKB are increased in the anaplastic thyroid carcinoma cell lines ARO, Cal 62, and KAT-4(PMID:15562011).

Notable Publications

Author	Pubmed ID	Journal	Application
Weiye Zhu	38142764	Exp Eye Res	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data