

For Research Use Only

Mhc Polyclonal antibody

Catalog Number: 10799-1-AP

15 Publications



Basic Information

Catalog Number:

10799-1-AP

Concentration:

900 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1254

GenBank Accession Number:

NM_078863

GeneID (NCBI):

35007

Full Name:

Myosin heavy chain

Calculated MW:

221 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

Cited Species:

mouse, rat

Background Information

The function of Mhc is muscle contraction. This antibody detects the Myosin heavy chain of drosophila.

Notable Publications

Author	Pubmed ID	Journal	Application
Cheng Chen	33007290	Eur J Pharmacol	WB
Ying Li	34548056	BMC Musculoskelet Disord	WB
Yuxiang Wang	35881904	J Cardiovasc Pharmacol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

