

For Research Use Only

# DDX5,p68 Polyclonal antibody

Catalog Number: 10804-1-AP

Featured Product

10 Publications



## Basic Information

**Catalog Number:**

10804-1-AP

**Size:**

750 µg/ml

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG1237

**GenBank Accession Number:**

BC016027

**GeneID (NCBI):**

1655

**UNIPROT ID:**

P17844

**Full Name:**

DEAD (Asp-Glu-Ala-Asp) box polypeptide 5

**Calculated MW:**

69 kDa

**Observed MW:**

65-68 kDa

**Purification Method:**

Antigen affinity purification

**Recommended Dilutions:**

WB 1:500-1:1000

IHC 1:20-1:200

## Applications

**Tested Applications:**

IHC, WB, ELISA

**Cited Applications:**

WB, IP, IF, IHC

**Species Specificity:**

human, mouse, rat

**Cited Species:**

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

**Positive Controls:**

WB : NIH/3T3 cells, A431 cells, HeLa cells

IHC : human breast cancer tissue,

## Background Information

DDX5 or p68, a member of DEAD/H BOX, which has a characteristic DEAD box, is a proliferation-associated nuclear antigen. DDX5 involves in the alternative regulation of pre-mRNA splicing, act as a RNA helicase to increase tau exon10 inclusion in a RBM4-dependent manner. It's also a transcription coactivator for estrogen receptor ESR1 and androgen receptor AR. Once synergized with DDX17 and SRA1 RNA, DDX5 activated MYOD1 transcription activity and involved in skeletal muscle differentiation. For the reason that DDX5 can be phosphorylated and dephosphorylated at several sites, at some cases the molecular weight of DDX5 is less than 68kD, approximately 65kd

## Notable Publications

Author	Pubmed ID	Journal	Application
Norihiro Okada	34082065	Gene	WB
Lingcai Zhao	35583334	Microbiol Spectr	WB,IF,IP
Ping Qiu	35672291	Cell Death Discov	WB,IF

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 10804-1-AP (DDX5, p68 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 10804-1-AP (DDX5, p68 Antibody) at dilution of 1:50 (under 40x lens).