

For Research Use Only

# Collagen Type IV Polyclonal antibody

Catalog Number: 10807-1-AP

7 Publications



## Basic Information

Catalog Number:

10807-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1118

GenBank Accession Number:

BC005305

GeneID (NCBI):

1288

UNIPROT ID:

Q14031

Full Name:

collagen, type IV, alpha 6

Calculated MW:

1691 aa, 164 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse

Cited Species:

human, mouse, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yuqi Li	36132982	J Immunol Res	IHC
Xiang-Hua Yu	36321394	Oral Dis	IF
G X Feng	28846110	Oncogene	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

