

For Research Use Only

Collagen Type IV Polyclonal antibody



Catalog Number: 10807-1-AP

5 Publications

Basic Information

Catalog Number:

10807-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1118

GenBank Accession Number:

BC005305

GeneID (NCBI):

1288

UNIPROT ID:

Q14031

Full Name:

collagen, type IV, alpha 6

Calculated MW:

1691 aa, 164 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, mouse

Cited Species:

human, rat, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuqi Li	36132982	J Immunol Res	IHC
Xiang-Hua Yu	36321394	Oral Dis	IF
G X Feng	28846110	Oncogene	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data