

For Research Use Only

PTHLH Polyclonal antibody

Catalog Number: 10817-1-AP

10 Publications



Basic Information

Catalog Number:

10817-1-AP

Size:

270 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1133

GenBank Accession Number:

BC005961

GeneID (NCBI):

5744

UNIPROT ID:

P12272

Full Name:

parathyroid hormone-like hormone

Calculated MW:

20 kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Moyuan Deng	34621738	Front Cell Dev Biol	WB
Qiusheng Shan	36238359	Oncol Lett	IHC
Kazumichi Kitayama	34713296	Int J Oncol	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

