

For Research Use Only

PTN Polyclonal antibody

Catalog Number: 10821-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

10821-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1151

GenBank Accession Number:

BC005916

GeneID (NCBI):

5764

UNIPROT ID:

P21246

Full Name:

pleiotrophin

Calculated MW:

18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhenyu Wang	29073605	Cell Physiol Biochem	WB
Zikai Zhou	34407403	Cell Rep	
Weiwen Wang	28447229	Angiogenesis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

