

For Research Use Only

DUSP9 Polyclonal antibody

Catalog Number: 10826-1-AP

6 Publications



Basic Information

Catalog Number:

10826-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1259

GenBank Accession Number:

BC060837

GeneID (NCBI):

1852

UNIPROT ID:

Q99956

Full Name:

dual specificity phosphatase 9

Calculated MW:

42 kDa

Observed MW:

42-44 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qiang Wei	31772940	J Diabetes Res	WB, IHC
Yuki Kanda	33249692	Genes Cells	WB
Lang Jiang	34239349	Int J Biol Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

