

For Research Use Only

PCDHGC3 Polyclonal antibody

Catalog Number: 10857-1-AP



Basic Information

Catalog Number:

10857-1-AP

Concentration:

247 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1243

GenBank Accession Number:

BC019299

GeneID (NCBI):

5098

UNIPROT ID:

Q9UN70

Full Name:

protocadherin gamma subfamily C, 3

Calculated MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

IHC : human breast cancer tissue,

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer using 10857-1-AP (PCDHGC3 antibody) at dilution of 1:100 (under 10x lens).