

For Research Use Only

LSM10 Polyclonal antibody

Catalog Number: 10874-1-AP



Basic Information

Catalog Number:

10874-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1027

GenBank Accession Number:

BC007623

GeneID (NCBI):

84967

UNIPROT ID:

Q969L4

Full Name:

LSM10, U7 small nuclear RNA

associated

Calculated MW:

14 kDa

Observed MW:

14 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human spleen tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human spleen tissue were subjected to SDS PAGE followed by western blot with 10874-1-AP (LSM10 antibody) at dilution of 1:200 incubated at room temperature for 1.5 hours.