

For Research Use Only

# CRK Polyclonal antibody

Catalog Number: 10973-1-AP



## Basic Information

Catalog Number:

10973-1-AP

Size:

330 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1420

GenBank Accession Number:

BC009837

GeneID (NCBI):

1398

UNIPROT ID:

P46108

Full Name:

v-crK sarcoma virus CT10 oncogene  
homolog (avian)

Calculated MW:

40 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate cancer using 10973-1-AP (CRK antibody) at dilution of 1:50 (under 10x lens).