

For Research Use Only

GRP94 Polyclonal antibody

Catalog Number: 10979-1-AP

Featured Product

7 Publications



Basic Information

Catalog Number:

10979-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1439

GenBank Accession Number:

BC009195

GeneID (NCBI):

7184

UNIPROT ID:

P14625

Full Name:

heat shock protein 90kDa beta (Grp94), member 1

Calculated MW:

96 kDa

Observed MW:

94 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

IF-P 1:50-1:500

IF/ICC 1:200-1:800

Applications

Tested Applications:

WB, IHC, IF/ICC, IF-P, FC (Intra), IP, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat, pig

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : A549 cells, HEK-293 cells, K-562 cells, MCF-7 cells, NIH/3T3 cells

IP : MCF-7 cells,

IHC : human colon tissue, human breast cancer tissue, human lung cancer tissue, mouse brain tissue

IF-P : mouse pancreas tissue, human appendicitis tissue

IF/ICC : HEK-293 cells,

Background Information

Glucose-regulated protein 94 (GRP94), also called HSP90B1 and GP96, is an endoplasmic reticulum (ER)-resident member of the heat shock protein 90 (HSP90) family (PMID: 33802964; 24658275). Under ER stress conditions, GRP94 accelerates its function as a molecular chaperone. Client proteins of GRP94 include toll-like receptors (TLRs), glycoprotein (GP) IX subunit, insulin-like growth factors (IGFs), proinsulin, and integrins (PMID: 7913987; 32781621; 22079671). As well as being an ER chaperone, GRP94 can also participate in Calcium regulation and is also a regulator of innate and adaptive immunity (PMID: 33802964; 11584270).

Notable Publications

Author	Pubmed ID	Journal	Application
Maojun Liu	28990064	Mol Med Rep	WB
Liang Zhang	33721634	Vet Microbiol	WB
Na Liu	35093118	J Neuroinflammation	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of GPR94 antibody (10979-1-AP, 1:1000) with si-control and si-GRP94 transfected HEK293 cells..




Various lysates were subjected to SDS PAGE followed by western blot with 10979-1-AP (GRP94 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human breast cancer using 10979-1-AP (HSP90B1 antibody) at dilution of 1:50 (under 10x lens).




Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using GRP94 antibody (10979-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).




IP result of anti-GRP94 (IP:10979-1-AP, 4ug; Detection:10979-1-AP 1:500) with MCF-7 cells lysate 4000ug.




Immunofluorescent analysis of (4% PFA) fixed mouse pancreas tissue using GRP94 antibody (10979-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



1×10^6 MCF-7 cells were intracellularly stained with 0.8 ug GRP94 Polyclonal antibody (10979-1-AP) and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.8 ug Rabbit IgG control Rabbit PolyAb (30000-O-AP) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 10979-1-AP (GRP94 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 10979-1-AP (GRP94 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).