

For Research Use Only

TSNAX Polyclonal antibody, PBS Only

Catalog Number:11040-1-PBS



Basic Information

Catalog Number:

11040-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1507

GenBank Accession Number:

BC011797

GeneID (NCBI):

7257

UNIPROT ID:

Q99598

Full Name:

translin-associated factor X

Calculated MW:

33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human gliomas using 11040-1-AP (TSNAX antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 11040-1-PBS in a different storage buffer formulation.



IP result of anti-TSNAX (IP:11040-1-AP, 4ug; Detection:11040-1-AP 1:300) with HeLa cells lysate 2400ug. This data was developed using the same antibody clone with 11040-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 11040-1-AP (TSNAX antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11040-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using TSNAX antibody (11040-1-AP) at dilution of 1:400 and CoralLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 11040-1-PBS in a different storage buffer formulation.