

For Research Use Only

LUZP4 Polyclonal antibody

Catalog Number: 11072-1-AP



Basic Information

Catalog Number:

11072-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1539

GenBank Accession Number:

BC128134

GeneID (NCBI):

51213

UNIPROT ID:

Q9P127

Full Name:

leucine zipper protein 4

Calculated MW:

40 kDa

Observed MW:

36 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:500

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MDA-MB-231 cells, SKOV-3 cells, U-87 MG cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11072-1-AP (LUZP4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.