

For Research Use Only

CLK2 Polyclonal antibody

Catalog Number: 11092-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

11092-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1549

GenBank Accession Number:

BC014067

GeneID (NCBI):

1196

UNIPROT ID:

P49760

Full Name:

CDC-like kinase 2

Calculated MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : A549 cells, Raji cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yongshuai Jin	26124922	Genes Cancer	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11092-1-AP (CLK2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



WB result of CLK2 antibody (11092-1-AP; 1:1200; incubated at room temperature for 1.5 hours) with sh-Control and sh-CLK2 transfected A549 cells.