

For Research Use Only

# Filamin A Polyclonal antibody

Catalog Number: 11140-1-AP

1 Publications



## Basic Information

Catalog Number:

11140-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1581

GenBank Accession Number:

BC014654

GeneID (NCBI):

2316

UNIPROT ID:

P21333

Full Name:

filamin A, alpha (actin binding protein 280)

Calculated MW:

2647 aa, 280 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

### Notable Publications

| Author       | Pubmed ID | Journal        | Application |
|--------------|-----------|----------------|-------------|
| Weiwei Liang | 38614124  | Endocr Connect | WB          |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

